

The enclosed vial may be analyzed for BOD, CBOD, COD and TOC.

Use a volumetric pipet to transfer 20 ml of the sample to a one (1) Liter volumetric flask and dilute to volume with laboratory grade water. Samples must be seeded for Biochemical Oxygen Demand and CBOD analysis. pH adjustment is required.

Analyte	Units	DMRQA
Biochemical Oxygen Demand	mg/L	9.00 - 158
Carbonaceous BOD	mg/L	8.00 - 136
Chemical Oxygen Demand	mg/L	15.0 - 256
Total Organic Carbon	mg/L	6.00 - 101

